

## **Self-heating compress**



- Ideal for home use
- □□ Only lumbar model available

It heats itself and is reusable. By simply tapping the inner pad, the compress reaches a temperature of approximately 50° for 30 minutes. Put in boiling water so that it regains its initial properties.

## 2 models avalable:

- Lumbar model. Dimensions : 36 x 13.5 cm. Adapted to the lumbar belt.
- Cervical model. Dimensions :  $42 \times 20$  cm. Ideal for the back of the neck. (This model is no longer available for sale).

## Variantes:

- Cervical
- Lumbar







- Lumbar - 835168

**Date d'impression :** 30/04/2025