

Self-heating compress



- Ideal for home use
- □□ Only lumbar model available

It warms up on its own and is reusable. Simply tap the internal pad to bring the compress to approximately 50°C for 30 minutes. Place in boiling water to restore its original properties.

2 models available:

- Lumbar model. Dimensions: $36 \times 13.5 \text{ cm}$. Suitable for lumbar belts.
- Cervical model. Dimensions: $42 \times 20 \text{ cm}$. Ideal for the neck. (This model is no longer available for sale)

Variantes:

- Cervical
- Lumbar







- Lumbar - 835168

Date d'impression : 02/11/2025