

## Giro standard wheelchair


















Ideal for crossing narrow door thresholds

The only chair that can also be used without the

The Giro chair has 2 types of wheels as standard. With its 24-inch wheels, it can be used like a classic chair. However, when you go home, for example, the helper can engage the small wheels underneath the wheelchair and remove the 24-inch wheels without having to lift the user.

By removing the 24-inch wheels, you can handle a chair with a reduced width (only 58 cm) and to pass door thresholds at home more easily. The small wheels are very stable and make the chair very manoeuvrable.

24 inch wheels with quick release.

24-inch wheels.

Removable armrests to facilitate lateral transfers. Height-adjustable, removable and retractable footrests. The seat and backrest covers are sewn together avoiding the gap between the seat and backrest. Anodised aluminium frame.

Seat dimensions: Width according to model x depth 41 cm x height 51 cm.

Backrest height 42 cm.

Handle height 91 cm.

Overall depth 107 cm.

Armrest height 20 cm.

Front wheels 20 cm.

Width of wheels 24" according to model: 60cm / 64cm /

67cm

Width small wheels according to model:  $54 \text{cm} \ / \ 58 \text{cm} \ / \ 61$ 

Folded width all models: 35,5 cm

Weight according to model: 15 kg / 15.4 kg / 16 kg.

Resistance: 140 kg.

Chassis guaranteed for 2 years.

3 different seat widths: 41, 45 or 48 cm (This size is no longer available for sale).

2 different colours to choose from : Pistachio or grey (Model grey in 45 cm is no longer available for sale). Model standard available with third party brakes.

## Variantes:

- 41 cm
- 45 cm (18 inches)
- 48 cm
- Grey
- Pistachio

**Date d'impression :** 02/08/2025