

**Surfa'safe Premium diffuse foam 750 ml**



Alcohol-free disinfectant foam for cleaning and disinfecting surfaces and non-submersible, non-invasive medical devices. Effective against microbes in just 2 minutes, dries quickly and leaves no residue. Ready to use.

Microbiological properties : Bacteria: EN 1040, EN 13727, EN 13697: 2 minutes. According to EN 13727: active against ERV, A. baumannii  $\beta$ LSE, E. coli OXA 48, E cloacae OXA 48, K. pneumoniae OXA 48, MRSA, Listeria monocytogenes, Salmonella enteritidis : 1 minute. According to EN 13697: active on ERV, A. baumannii  $\beta$ LSE, E. coli OXA 48, E cloacae OXA 48, MRSA, Listeria monocytogenes, Salmonella enteritidis: 2 mins. Yeasts : EN 1275, EN 13697: 5 mins. EN 13624: 2 mins. Mould: EN 1275, EN 13624, EN 13697 : 20 minutes. According to EN 13624, EN 13697 : Active against Aspergillus fumigatus and Tricophyton mentagrophytes: 15 minutes. Viruses: EN 14476 : Adenovirus, Murine Norovirus : 30 minutes. According to EN 14476 : active against HIV-1, PRV (HBV model virus), BVDV (HCV model virus), Vaccinia virus, HSV, Rotavirus: 1 min. According to EN 14476 : active against Polyomavirus SV 40, Coronavirus: 5 min.

Dangerous - follow the precautions for use. Use biocides with caution. Class IIa medical device. Before use, read the label and product information.

750 ml bottle with foam dispenser.

**Date d'impression : 27/06/2026**