

## Anios'Clean Excel D



- **Sanitises instruments**

Cleaning and disinfects medical devices : surgical and medical instruments, heat-sensitive instruments and endoscopy equipment, prior to sterilisation.  
Product to be diluted : 25 ml per 5 litres.

Microbiological properties : Bacteria : EN 13727, EN 14561 : 5 mins - 0.5%. EN 13727 on BMR : 15 min - 0.5%. Mycobacteria : EN 14348, EN 14563 (M. terrae) : 60 min - 2%. Yeasts : EN 13624, EN 14562 : 15 min - 0.5%. Viruses : PRV (HBV model virus), BVDV (HCV model virus), Herpesvirus, Vaccinia virus, HIV-1 : 5 minutes - 0.5%.

Dangerous - follow the precautions for use. Class IIb medical device. - Before use, read the label and product information.

Available in 1 or 5 litre containers.

**Variantes :**

- 1 L
- 5000 ml