

## Steranios 2%



- **High-level disinfection of medical devices, surgical and medical equipment, endoscopy equipment and heat-sensitive equipment**

High-level disinfection of medical devices, surgical and medical equipment, endoscopy equipment and heat-sensitive equipment.

Microbiological properties : Bacteria: EN 1040, EN 13727, EN 14561: 5 minutes. Active against *Helicobacter pylori*: 10 minutes. Mycobacteria: *Mycobacterium tuberculosis* (B.K.): 5 minutes. EN 14348, EN 14563 (*M. terrae*): 10 minutes. Active against *Mycobacterium massiliense* : 10 minutes. Yeasts/moulds: EN 1275. EN 13624, EN 14562: 10 minutes. Viruses: EN 14476: Poliovirus, Adenovirus, Murine Norovirus: 10 minutes. Active against HIV-1, HBV, BVDV (HBV model virus) and HSV: 5 minutes. Spores: NF T 72-230/231 *Bacillus subtilis*, *Bacillus cereus*, *Clostridium difficile*.

Dangerous - follow the precautions for use. Class IIb medical device. - Before use, read the label and product information.

Ready to use.

Available in 5-litre containers.

**Date d'impression : 27/04/2026**