

Fifth wheel Blumil GO



- 875030

- **Unique, patented fastening mechanism**
- **Quality of manufacture and safety in use**
- **Optional extra battery available**

The Blumil GO motorised system transforms your manual wheelchair into a power-assisted chair. The unique fastening system does not require any special adaptation. At the touch of a button, the connecting bar is attached to the axle of your wheelchair and electrically tilts the seat backwards. It is then ready to use. Simply position the footrests of your wheelchair on the floor support of the Blumil GO. You can also switch it off in a matter of seconds. Equipped with a 20 cm wheel for use in the city and a USB charging socket. Electronic display. Double LED lighting. 300W motor. Battery LG 7,8 Ah 280 Wh. Total weight with battery 13 kg. Max. speed 10 km/h. Autonomy 30 kms.

Optional extra battery available

Product compatible with active chairs only