

MOOVY ELECTRIC TOILET / SHOWER AND TRANSFER CHAIR

• The electric lift version of the Moovy chair

The electrically operated version of the Moovy Transfer Toilet / Shower Chair.

Like the standard Moovy, it allows caregivers to transfer users easily, effortlessly and without risk of pain. Transfer users from bed to shower or toilet.

Simply position the Moovy chair at the height of the user's buttocks, adjusting the height using the electric hydraulic actuator. The up/down control is located directly on the Moovy Chair (no remote control). Then open the seat, close it around the user and secure it with the double safety mechanism (belt clip and buckle).

It has four wheels with brakes, a flexible safety belt and a soft foam backrest and armrests.

It comes with a comfort cushion and a removable potty for independent use.

It also has foldable handles for easier access to the patient.

It can be used in the bathroom or shower provided that the battery is removed before use. Please consult the instructions and follow the precautions for use and operating instructions for more information.

Overall dimensions : width 62.5 x depth 78 cm (floor) or 82 cm (top).

Seat dimensions : width 47 x depth 43 x height 44.5 to 64 cm.

Backrest height : 40.5 cm.

Width between armrests : 45.5 cm.

Weight : 35 kg.

Maximum user weight : 150 kg

