

Stairlift chair



- **Go up and down the stairs in complete safety**

Simply place the user on the comfortable seat and secure them with the 3-point safety belt, footplate and armrests. By simply activating the handles, the stairlift starts moving smoothly and safely. The whole unit can be dismantled into 3 parts, making it easy to transport.

Self-sufficient, it runs on a rechargeable lithium 48V 10aH battery and a 48V 400W motor. Range: between 1,450 and 1,800 steps.

Choice of 12 or 17 steps per minute. Detachable motor for use on flat surfaces while pushing the wheelchair. Quiet, less than 70 dB in operation.

Dimensions: width 61.5 x depth 91.8 x height 151 to 179 cm (adjustable handles).

Weight: 37 kg.

Seat dimensions: width 47 x depth 48 x height 52 cm.

Resistance 140 kg