

## Norco glove



- Wide right hand - 116007
- Right hand medium - 116005
- Right hand small - 116003
- Large left hand - 116006
- Left hand medium - 116004  
(On restock - Available from 09/05/2024)
- Left hand small - 116002

- **Compression glove that compresses the tissues in order to control the development of oedemas.**
- **Models for left or right hand with 3 sizes per hand**

Our lycra gloves gently compress the tissues to control the development of oedema. They help to avoid certain joint pains in arthritic patients. The fabric used allows air to circulate and does not cause heating.

Compression models for left or right hand with 3 sizes per model.

Measure the circumference of your hand at the level of the phalanges at the base of the fingers to choose your size: small (18 to 20 cm), medium (19 to 22 cm) and large (20 to 25 cm). Compression 15 to 25 mm Hg

Care instructions: Machine wash cold (delicate cycle). Dry flat in the open air. If there is too much thread, do not cut it to avoid fraying the seam.

### **Variantes :**

- Wide right hand
- Right hand medium
- Right hand small
- Large left hand
- Left hand medium
- Left hand small