


## **Capiluve shampoo basin for bedridden patients**



ė.



- 823059

Shampoo basin for easy and effortless shampooing of bedridden patients

The Capiluve makes it very easy to wash the hair of bedridden patients. Simply slide it under the patient's shoulders so that it stays firmly in place, and position the neck in the basin. The wide basin allows for freedom of movement. Water does not stagnate, but flows into the Capiluve's internal reservoir, which can hold up to 8 litres. When the shampoo is finished, the entire unit can easily be transported without any water leakage and emptied using the drain plug. Effortless shampooing made easy.

Dimensions :  $60 \times 60 \times 16 \text{ cm}$ .

Weight: 2.7 kg.

Material: Plastic polyethylene

**Date d'impression :** 19/08/2025