

Capiluve shampoo basin for bedridden people



- 823059

- **Container for easy and effortless shampoo of the bedridden person**

Our shampoo basin Capiluve is designed to shampoo the head of bedridden people very easily. Simply slide it under the patient's shoulders so that it holds firmly, and position the back of the neck in the basin. The wide pelvis gives freedom of movement. The water does not stagnate. It flows into the internal tank of the capiluve, which can hold up to 8 litres. When the shampoo is finished, the whole unit can be easily transported, without any water leakage, to be emptied via the drain plug. An easy shampoo without unnecessary effort.

Dimensions: 60 x 60 x 16 cm.

Weight 2.7 kg.

Material : Polyethylene plastic

Date d'impression : 03/05/2024