

Victor transfer rollator



- **2 in 1 : rollator and transfer chair**
- **With comfortable armrests**

The Victor rollator can easily be converted into a transfer chair by unfolding the footrest and changing the direction of the backrest. Ideal for walking safely and being brought back to a seated position when tired. Equipped with a padded seat, handles, and a foam backrest. Simple braking system + parking brake.

Large 20 cm diameter wheels.

Overall dimensions : width 60 x depth 65 x height 78 to 89 cm in 6 positions.


Seat dimensions : width 35 x depth 32 x height 52 cm.

Width between armrests : 45.5 cm.

Weight : 8.8 kg.

Color : burgundy.

Maximum user weight : 110 kg.



Étapes pour inverser le dossier



Étapes pour monter les accoudoirs



Étapes réglage hauteur des poignées



Date d'impression : 19/12/2025