

Rollator transfer chair Victor



- **2 in 1: rollator and transfer chair**
- **With comfort armrests**

Our Victor rollator easily converts to a transfer chair, simply unfold the footrest and change the direction of the backrest. Ideal for walking in safety and to be brought back seated in case of fatigue. Equipped with a padded seat, handles and foam backrest. Simple braking system + parking.

Large wheels with a diameter of 20 cm.
Overall dimensions : width 60 x depth 65 x height 78 to 89 cm in 6 positions.

Seat dimensions: width 35 x depth 32 x height 52 cm.

Width between armrests : 45,5 cm.

Weight 8.8 kg.

Color burgundy.

Maximum weight 110 kg.



Étapes pour inverser le dossier



Étapes pour monter les accoudoirs



Étapes réglage hauteur des poignées



Date d'impression : 30/04/2025