

## Heated lumbar belt



- Self-heating - 835180

- Thermo-cryo - 835181

- **3 models available**

Very easy to use, simply attach the lumbar belt to the waist to benefit from the heat on the lumbar region. This belt can also be used wrapped around a thigh or an arm.  
Length: 120 cm.

3 models available:

- Hot cold model with microwave heating  
(This model is no longer available for sale).
- Auto model with automatic heating
- Thermo-cryo model with bowl heating.

**Variantes :**

- Self-heating
- Hot Cold
- Thermo-cryo