

Malte suspended toilet armrests



- **No obstacles on the floor**
- **Sturdy and stable**

The Malte toilet armrests can be installed on most toilet bowls. They make it much easier to sit down and stand up, and provide stability and support. They can be installed under the existing toilet seat or under a toilet seat raiser. They are adjustable in width with two positions to suit all body types.

Integrated toilet roll holder.

Steel frame and foam handles.

Dimensions : width 48 x depth 24 cm.

Height of armrests above the toilet bowl : 26 to 38 cm in 6 positions.

Centre-to-centre distance : 15 cm.

Maximum user weight : 180 kg.

VIDÉO DISPONIBLE
JUSTE EN DESSOUS !



- 811162