

Pic Enema bulb



- 35 ml - 802067
- 70 ml - 802068
- 125 ml - 802363
- 365 ml - 802364

- **With soft cannula : A gentle and comfortable cleansing**
- **With a rigid cannula: Ideal for simple bowel care**

[Pic Enema bulb](#) with rigid or soft cannula depending on the model. Suitable for auricular or intestinal cleansing depending on the model :

- Flexible: Rubber bulb with a flexible tube. Suitable for colonic irrigation and ear washing (small size only).
- Rigid: Bulb for colonic irrigation, made from a sturdy, soft and non-toxic material.

Available in several models:

- 35 ml. With soft cannula.
- 70 ml. With soft cannula.
- 125 ml. With soft cannula.
- 365 ml. With rigid cannula (9 cm).

Enema spray: PVC. Non-sterile.

Before each use, wash the device with soap and water and rinse thoroughly. After use, wash with soap and water and rinse thoroughly. Allow all components to dry before storing. Before each re-use, repeat the operations specified from the beginning.

Variantes :

- 35 ml
- 70 ml
- 125 ml
- 365 ml