

Hebix lifting aid chair



VIDÉO DISPONIBLE
JUSTE EN DESSOUS !



- Dismountable chair - 821040
- Transport bag - 821106
- Chest strap - 821109
- Folding chair - 821041

- **A simple and innovative solution to gently lift people off the floor**
- **Ideal to relieve the caregivers**

The Hebix can be opened in a few simple steps and placed over the person lying on the floor. The key factor here is that the force required is only about a quarter of the body weight of the person being lifted : for someone weighing 80 kg, this is only about 20 kg. By simply lifting the Hebix using the handles, it is possible to gently and painlessly raise the person lying down into a sitting position. They can remain in this position, stand up or be moved – it couldn't be easier!
2 models available :

- Model that can be dismantled into 3 parts. Weight : 7 kg.
- Foldable model is more suitable for emergency services. Weight : 10.5 kg.

Dimensions : Length 75 x width 45 cm.
Maximum user weight : 100 kg.

Optional carry bag available for the foldable model only.
Optional chest strap for added safety during lifting.

Variantes :

- Dismountable chair
- Transport bag
- Chest strap
- Folding chair