

Modulo 2-wheel rollator



- Removable comfortable seat
- 2 zipped pockets under the seat
- Height-adjustable seat

The Modulo 2-wheel rollator can be used without the inconvenience of a seat in the center of the frame, simply by lifting the seat and tilting it forward. When the user wishes to sit down, it is simple to replace it on the supports. The height of the handles can be adjusted using a lever and indexing system. Very sturdy aluminum frame.

Overall dimensions : width 59 x depth 55 x height 83 to 97 $\,$

cm in 5 positions.

Seat dimensions : width 39 x depth 34 x height 50 / 53 / 56

cm.

Weight: 4.8 kg.

Maximum user weight: 136 kg



























- 826038

Date d'impression: 05/11/2025