

GelPro hot-cold compresses - hand / wrist



- **Revolutionary stretch material !**
- **Enables efficient cold or hot treatment**

The new GelPro range offers complete thermotherapy treatment. Place the compress in the freezer to relieve minor injuries and swelling. Place the compress in the microwave to soften and relieve muscle or organ pain. The compresses slip onto the area to be treated and fit snugly because the material used is elastic, allowing it to adapt to most body types.

Colour : blue.

Model slips over the palm to the wrist, with separate thumb.
Dimensions without extension : 17 x 10.5 cm.

Compresses available in 5 other models.



- 835187
(On restock - Available from 19/11/2025)