

## GelPro hot-cold compresses - shoulder



- Revolutionary stretch material!
- Enables efficient cold or hot treatment

The new GelPro range offers complete thermotherapy treatment. Place the compress in the freezer to relieve minor injuries and swelling. Place the compress in the microwave to soften and relieve muscle or organ pain. The compresses slip onto the area to be treated and fit snugly because the material used is elastic, allowing it to adapt to most body types.

Colour : blue.

Half belt design with tube to slip over the arm to cover the entire shoulder. The elastic velcro strap fastens under the opposite arm.

Dimensions without extension : 29 x 35 cm.

Strap length without extension : 76 cm.

Compresses available in 5 other models.



- 835192

(On restock - Available from 10/02/2026)

**Date d'impression :** 02/02/2026