

## Complete chest support kit



- Size 1 - For rigid backrests - 833031.T1
- Size 1 - For flexible backrests - 833030.T1
- Size 2 - For rigid backrests - 833031.T2
- Size 2 - For flexible backrests - 833030.T2
- Size 3 - For rigid backrests - 833031.T3
- Size 3 - For flexible backrests - 833030.T3

- **2 kits to choose from.**
- **Consists of a canvas and 2 side supports to be positioned as required**
- **Supports can be repositioned quickly**

Physipro chest supports provide lateral support for the upper body, from the hip bone to the armpit. Highly efficient in cases of hypotonia, they can be adapted to all existing chairs (wheelchairs or geriatric chairs).

### Quick and easy to install :

The fabric is secured to the backrest using three buckle straps and four straps at the corners. The supports can then be positioned at the desired height and angle. They can be repositioned as required.

### Adaptable to your existing chairs :

- On flexible backs (wheelchairs) : choose the anti-hammock effect kit because the angle of the supports is open to compensate for the distortion of the fabric. Consists of a fabric and 2 side supports. Choice of sizes.

- On rigid backs (geriatric chairs) : choose the kit with a 90° angle of the supports.

Consists of a fabric cover and two side supports. Choice of sizes.

### Available sizes :

- Size 1: Width 35 cm x height 40 cm, suitable for chairs 35 to 40 cm wide
- Size 2: Width 40 cm x height 40 cm, suitable for chairs 40 to 45 cm wide
- Size 3: Width 45 cm x height 40 cm, suitable for chairs 45 to 50 cm wide

Standard depth 16.75 or 13.75 cm on request.

### Variantes :

- Size 1
- Size 2
- Size 3
- For rigid backrests
- For flexible backrests